INSTRUCTION MANUAL
BEDIENUNGSANLEITUNG
MANUEL D'UTILISATION
MANUALE D'ISTRUZIONI

Rotator-Mixer RM Multi-1

N2400-5000

www.starlabgroup.com

CE MARKING

The CE marking affixed to the equipment indicates that the

STARLAB Rotator-Mixer RM Multi-1
Cat. No. N2400-5000

meets the requirements of the following directives:

ELECTROMAGNETIC COMPATIBILITY

Council Directive 89/336/EEC

Applied Standards: **EN 61326**

Electrical equipment for measurement,
control and laboratory use, EMC Requirements

COUNCIL DIRECTIVE 73/23/EEC

**Electrical Equipment Designed
for the Use Within Certain Voltage Limits**

Applied Standards: **EN 61010**

Electrical equipment for measurement,
control and laboratory use, EMC Requirements

STARLAB GmbH

| | |
|----------------|-----------|
| English | 3 |
| German | 11 |
| French | 19 |
| Italian | 26 |

TABLE OF CONTENTS

| | |
|--------------------------------|----|
| Safety Precautions..... | 3 |
| 1. General Information..... | 4 |
| 2. Getting Started | 5 |
| 3. Operation | 6 |
| 4. Program setting | 7 |
| 5. Specifications | 9 |
| 6. Technical Maintenance | 9 |
| 7. Warranty | 10 |

SAFETY PRECAUTIONS

Please read these operating instructions fully before use.


- Use only as specified within these operating instructions.
- Do not use power supply units other than those recommended by the manufacturer.
- Connect only to a power supply with a voltage that corresponds to that stated on the serial number label.
- Ensure that the mains switch and isolating device (power supply connector) are easily accessible during operation.
- To turn off the device, disconnect the power supply unit from the power socket.
- Do not operate the unit near aggressive or explosive chemical mixtures.
- When loading, use an even number of tubes arranged symmetrically to the rotation axis, to ensure balance during operation.
- Check with the manufacturer before using any cleaning or decontamination method that is not recommended to ensure it will not damage the equipment.
- Do not operate the unit outside the laboratory premises.
- Disconnect at the power supply socket before moving.
- If liquid is spilt inside the unit, disconnect it from the power supply. The unit should then be checked by a competent person.

1. GENERAL INFORMATION


The STARLAB programmable Rotator–Mixer RM Multi-1 provides orbital rotational motion, reciprocal motion, and vibro motion of the platform. The microprocessor control allows the setting of programs for both separate mixing motions and cycles of consecutive, different motion types.

SETTING OPTIONS:


1. Speed (1–100 rpm) and time of ordinary **ROTATIONAL MOTION** of the platform for a selected time period (0 to 250 seconds) or non-stop


2. Segment of **RECIPROCAL MOTION** when the direction of platform rotational motion is changing in turns within the limits of the set segment (turning angle 15–90°) for a selected time period (0 to 250 seconds) or non-stop

3.1 Segment and time of **VIBRO MOTION** of the platform (turning angle 1–5° for a time period 1–5 seconds) run on the borders of reciprocal motion segment. Available only when the reciprocal motion is ON.

3.2 PAUSE duration for temporary platform motion stoppage (1–5 seconds) when the vibro motion is off (the turning angle of vibro motion is set to zero) run on the borders of reciprocal motion segment. Available only when the reciprocal motion is ON.

4. Working period from 1 min to 24 hours, or non-stop.

The STARLAB Rotator–Mixer RM Multi-1 offers a user-friendly interface, which provides options not only for changing the program during the operation, but also for simultaneous control over different steps of mixing protocol realisation.

The STARLAB Rotator–Mixer RM Multi-1 will undoubtedly provide increased methodical means to researchers working in the field of modern molecular and cell biology and the developing biodiagnostics technology based on the use of magnetic particles, for which unexpected and disturbing hydrodynamic shifts of the reactants are essential.

The STARLAB Rotator–Mixer RM Multi-1 is designed for mixing applications including cell suspensions, incubations, extraction procedures, gel washing and precipitations.

The unit is applicable in all areas of laboratory research in biotechnology, microbiology, chemistry, and medicine.

2. GETTING STARTED

Unpacking:

Remove packing materials carefully and retain for future shipment or storage of the unit.

Check you have everything:

| | |
|----------------------------------|----|
| STARLAB Rotator–Mixer RM Multi-1 | x1 |
| Standard Platform | x1 |
| External Power Supply | x1 |
| Operating Manual | x1 |

An alternative optional platform is also available separately:

N2400-5101 for 4 x 50ml tubes and 12 x 1.5-15ml tubes.

Installation:

- Place the device in the working area on a horizontal, even working surface.
- Plug the external power supply unit into the 12V socket at the back of the Rotator.

Platform replacement:

Unscrew the two fixing screws on the platform. Replace the platform, and secure the new platform with the two screws.

3. OPERATION

Connect the external power supply to the mains.

Place tubes on the platform: Microcentrifuge tubes should be pushed into the rubber holder up to the lid (tube's maximum diameter); for vacutainers and tubes with caps, these should be pushed into the holder up to half the tube's length.

Set the appropriate program and operation time (see 4. PROGRAM SETTINGS).

Press the RUN/STOP key (1 / Fig.1) to start the program.

The platform motion begins, and the display shows the corresponding indication RUN (2 / Fig.1) and the changing time values.

If the operation time is not set and the timer indicator (3 / Fig.1) shows 0:00, pressing the RUN/STOP key will activate continuous mode until the Run/Stop key is pressed again.

If the operation time is set, the rotator will stop when the set time interval expires (flashing STOP on the display) and it will give a sound signal to indicate the end of operation (press RUN/STOP key to stop the signal).

To repeat the set program press RUN/STOP key.

If necessary, the Rotator can be stopped any time during operation before the set time is expired by pressing the RUN/STOP key. In this case, the platform motion stops when the platform achieves horizontal position. Pressing the RUN/STOP key again will start the program from the beginning (countdown timer will be restarted).

NOTE: A step motor is used in the Rotator. It is possible to stop the platform by hand for a moment without damaging the mechanical nodes of the device. If the platform is stopped by hand during operation, the program does not stop and the platform motion is automatically resumed after the platform is released.

Unplug the external power supply from the mains outlet to turn off the unit.


Fig. 1

4. PROGRAM SETTINGS

To change the parameter, press the SELECT key (4 / Fig. 1) (the active parameter will flash in the display), then use the UP and DOWN keys (5 / Fig. 1) to select the required parameter (Note: if the key is pressed and held for more that 2 seconds the numbers change quickly).

Program setting can be changed during operation. The updated/last setting will be automatically stored on the instrument's memory.

The countdown timer is used to control the operation time. The timer can be set for a time period of 1 minute to 24 hours (in increments of 1 minute).

The examples below show separate motion types and their available combinations in cycles. The data on the right-hand side show the possible parameter values for each type. Note that the speed range for the STARLAB Rotator is 1–100 RPM.

MOTION/ROTATION CYCLES:

Orbital

Set the speed of orbital rotation (1–100 RPM), the time of orbital rotation (1–250 secs), and the time for reciprocal motion to zero (off).

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 15-90° | 0-5° |
| 1-250 | OFF | OFF/1-5 |
| Time (sec.) | | |

Orbital + Reciprical

Set the speed (1-100 RPM) and time (1–250 secs) of orbital rotation.
Set the turning angle (15–19°) and time (1–250 secs) for reciprocal motion.
Switch off the vibro motion by setting the time of vibro motion to zero (off).

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 15-90° | 1-5° |
| 1-250 | 1-250 | OFF |
| Time (sec.) | | |

Orbital + Reciprical + Vibro

Set the speed (1-100 RPM) and time (1–250 secs) of orbital rotation.
Set the turning angle (15–19°) and time (1–250 secs) for reciprocal motion.
Set the turning angle (1–5°) and time (1–5 secs) for vibro motion

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 15-90° | 1-5° |
| 1-250 | 1-250 | 1-5 |
| Time (sec.) | | |

Note: if the set time of reciprocal motion is shorter than or equal to the time set for vibro motion, the reciprocal motion will be omitted resulting a in **Orbital + Vibro** cycle.

Orbital + Reciprocal + Pause

Set the speed (1–100 RPM) and time (1–250 secs) for orbital rotation.

Set the turning angle (15–19°) and time (1–250 secs) for reciprocal motion.

Set the turning angle of vibro to zero.

Set the time for Vibro/Pause mode (1–5 secs)

- this is the duration of the pause.

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 15-90° | 0° |
| 1-250 | 1-250 | 1-5 |
| Time (sec.) | | |

Note: if the set time of Reciprocal motion is shorter than or equal to the set time for Vibro/pause mode, the Reciprocal motion will be omitted resulting in a **Orbital + Pause** cycle.

Reciprocal Rotation

Set the speed (1–100 RPM) for orbital rotation.

Set the time for orbital rotation to zero (off).

Set the turning angle (15–19°) and time (1–250 secs) for reciprocal motion.

Set the time for Vibro to zero (off).

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 15-90° | 0-5° |
| OFF | 1-250 | OFF |
| Time (sec.) | | |

Reciprocal + Pause

Set the speed (1–100 RPM) for orbital rotation.

Set the time for orbital rotation to zero (off).

Set the turning angle (15–19°) and time (1–250 secs) for reciprocal motion.

Set the time for vibro motion (1–5 secs)

-this is the duration of the pause.

Set the angle of Vibro motion to zero (off).

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 15-90° | 0-5° |
| OFF | 1-250 | OFF |
| Time (sec.) | | |

Vibro + Reciprocal Rotation

Set the speed (1–100 RPM) for orbital rotation.

Set the time for orbital rotation to zero (off).

Set the turning angle (15–19°) and time (1–250 secs) for reciprocal motion.

Set the angle (1–5°) and time (1–5 secs) for Vibro motion.

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 15-90° | 0-5° |
| OFF | 1-250 | 1-5 |
| Time (sec.) | | |

Note: normally the rotator performs soft vibration (Vibro motion) but a mode for **hard vibration** is also available. To set hard vibration, set the turning angle of Reciprocal motion at 90° and the angle of Vibro motion to 1° (hard vibro).

| Orbital (rpm) | Reciprocal (deg.) | Vibro / pause |
|------------------|----------------------|------------------|
| 1-100 | 90° | 1° |
| OFF/1-250 | 1-250 | 1-5 |
| Time (sec.) | | |

5. SPECIFICATIONS

The unit is designed for operation in closed laboratory rooms at ambient temperature from +5°C to +40°C with a maximum relative humidity of 80% for temperatures up to 31°C, decreasing linearly to 50% relative humidity at 40°C.

| | |
|--|---|
| Speed Range | 1–100 RPM |
| Reciprocal Rotation Mode Turning angle Timer | 15–90° (increment 15°) 0–250 seconds |
| Vibro/Pause Mode Turning angle Timer Pause | 15–90° (increment 15°) 0–5 seconds 0–5 seconds |
| General Working Timer | 1 min to 24hrs (increment 1 minute), or non stop |
| External Power Supply | DC 12V, 500 mA |
| Standard Platform (supplied) | Holds 22 tubes up to 15mm diameter (1.5ml, 2.0ml, 15ml tubes) |
| Optional Platform Cat. No. N2400-5101 | Holds 4 tubes up to 30mm diameter (50ml tubes), and 12 tubes up to 15mm diameter (1.5ml, 2.0ml, 15ml tubes) |
| Dimensions (WxDxH) | 310 x 175 x 125mm |
| Maximum Load | 0.4 kg |
| Weight (excl. power supply) | 1.8 kg |

6. TECHNICAL

Maintenance and repairs should only be carried out by STARLAB Service personnel.

Standard ethanol (75%) or any other cleaning liquids which do not contain the strong organic solvents (such as acetone, toluene or benzene) can be used to clean and disinfect the Rotator.

7. WARRANTY

The STARLAB guarantee warrants the compliance of the Rotator–Mixer RM Multi-1 within the requirements of the Specifications, provided that the customer follows the operation, storage and transportation requirements

The Rotator is guaranteed for a period of 2 years from date of purchase.

If any manufacturing defects are discovered or you are not satisfied with the equipment, please contact your local Starlab office or distributor.

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STARLAB (UK), Ltd
STARLAB Srl
STARLAB India Ltd